


10. 12. 18

- 1) We turn on the laser and place glass point it at the glass with a randomly chosen angle. We place a Polarizer in its path and first deduce roughly the calibration by minimizing the intensity ~~that goes~~. (Around Brewster ~~angle~~ angle the polarisation is mainly S-Pol. \rightarrow Minimum corresponds to P-Pol)

$$\rho_1 \sim 94^\circ$$

$$\rho_2 \sim 120^\circ$$

pic 1: finding brewster



Use these settings to find brewster angle by minimizing intensity by changing angle of incident of the glass.


We find it around $54^\circ \approx \Theta_{\text{Brewster}}$

The Polarizer just lets light wave with a certain linear polarization through, meaning ~~by splitting the~~ the intensity will be reduced according to the amount of polarization. Beamsplitter can split the beam by just ~~as well as act like~~ halving the Intensity or splitting it according to

Intensity

- 2) If the second glass plate is installed with the x,y components of its normal being the same as the direction of incident and the angle between these is Θ_{Brewster} , the reflection disappears

pic 2:
Brewster
angle



no reflection cause the S-Polar of the first plate is P-Polar for the second \rightarrow Transmitted without reflection

- 3) To find the calibration we use same method as 1 just without changing glass plate angle, but instead polarization

ρ_1 ~~100,5°~~ 101° for $\leftarrow \rightarrow$ and 281° (\leftrightarrow) (P-Polari)


ρ_2 120° for $\leftarrow \rightarrow$ and 300° (\leftrightarrow) (P-Polari)

For the PBS you can compare its effect for PS-Polarized laser

and the "raw" laser to figure out which one does what.

→ one is 50/50 split
one is S/P-Polar split

The PBS doesn't change the Polarization.




Pic 3.
PBS Test
for Polar

11.12.18

ii) Check 3)

5)




Use this setup to calibrate $\lambda/2$
by minimizing Laser after S-Polarizer

$90^\circ, 180^\circ, 270^\circ, 360^\circ$ (a)

$104^\circ, 194^\circ, 284^\circ, 14^\circ$ (b)

→ If you keep rotating the $\lambda/2$ in the same direction ($360^\circ +$)
you move the calibration around $2-3^\circ$

6)



(a) $20^\circ, 200^\circ$ (b)
 $140^\circ, 320^\circ$ (a)

→ only 180° "repeatable" cause $\lambda/4$ -plate is "sensitive" to
direction of circular polarized light

7) Fig 3 ~~Fig 3~~ is same as pic 5

When u turn $2/4$ ~~at~~ 90° you reverse the direction
of polarized light and the second $2/4$ will turn

it into ~~go~~ turned initial linearized light

$$\begin{pmatrix} F_x \\ F_y \end{pmatrix} \rightarrow \begin{pmatrix} F_x + \frac{\pi}{2} \\ F_y \end{pmatrix}$$

$$L: \begin{pmatrix} F_x \\ F_y \end{pmatrix} = \begin{pmatrix} \cos \theta \\ \sin \theta \end{pmatrix}$$

$$\begin{pmatrix} -\sin \theta \\ \cos \theta \end{pmatrix}$$

$$\begin{array}{c} F \\ \downarrow \\ S \\ \downarrow \\ L \end{array}$$

$$\begin{pmatrix} -\cos \theta \\ \sin \theta \end{pmatrix}$$

Fig 4 and 3

The difference between



and




is none because only the relative position between light directions and Waveplate Fast/slow-Axis matters.

~~And if you turn a $\lambda/4$ -waveplate around.~~



Difference only shows when we swap the direction which the $\lambda/4$ -plate is facing. Here we swap the second $\lambda/4$ waveplate and deduce via a polarizer if the Polarization changed or not.



We put the for Fig 3 we put the first polarizer on $P_1 \cancel{235^\circ}$, $\lambda/4$ (a) 90° (c) 140° $P_2 240^\circ \Rightarrow$ No projection



We flip(a)



$P_2 300^\circ \Rightarrow$ No projection
Light is \perp polarized


For Fig 4

$P_1 235^\circ$ $\lambda/4$ (b) 65° $\lambda/4$ (a) 185° $P_2 255^\circ \Rightarrow$ No projection
Light is \perp polarized

Flip (a)



No projection
 \Rightarrow Light is still \perp polarized



So we set the 1. Polarizer to 180° to get an S-Pol.

~~But the second Polarizer is set to no "minimum" \rightarrow No linear Polarisation~~

We insert a $\lambda/4$ -Plate after mirror to test for circular Polarisation \rightarrow ? Polarizer sees minimum for which the 2. Poli can still find (minimum at $90^\circ \rightarrow$ Light is \perp)

Same for P-Poli. (we change it between 1. Poli and Mirror with $\lambda/4$ -Plate)

g)

We still get a linear Polarization even after the mirror between first and second $\lambda/4$ but the linear polarization is turned for roughly 100° .

After turning the mirror, by increasing the angle of reflection, we can see that the polarization is turning as well. One can see that Θ_r and Θ_p are



behaving similar. By increasing the Θ_r for roughly 40° the Polarisation (Θ_p) also turns 40° .

10)


The mirror should not have an influence on the polarization so that we look at the same experiment as in Fig 3 so if we expect it to be a linear Polarisation turned 90° in respect to its starting polarization.

We test it by having the mirror at almost 0° and placing $\lambda/4$ between the first Polarizer and the mirror.

The first Polarizer will give us \perp -polarized light but the reflected and light beam which passed through $\lambda/4$ disappears

Polarizer.

11) With the half $\lambda/2$ -Plate we expect the Polarization to ~~not~~ not turn since the $\lambda/2$ Plate will always turn if it's placed with 45° between its S/F-axis and PFS-Polarisator.




12)

~~The isolater tells~~

We check with the ceiling lights, one

One direction appears green, the other appears orange

The laser only is let through the isolater when it would show orange for white light.



Faraday-Effekt

13) It's the same

12.12.18

EOM

5.2.0.1

1/2-factor Power Supply/High Voltage

Wave Generator [V]	High Voltage [kV]
0,166	0,06
0,592	0,19
1,036	0,32
1,712	0,52
2,350	0,71
3,127	0,95
3,632	1,09
4,383	1,32
5,174	1,56
5,948	1,79

Label 1: Amplifier :


This equals roughly 300 as a factor as shown in Fig

→ Never exceed 5.5V to

below 1.8 kV


5.2.2

- 1) The optical axes are at $-90^\circ, 0, +90^\circ$ of the Pockel cell
- 2-3) Done \rightarrow Parallel beams, below 3V
- 4) $+90^\circ$ tek2.csv




5.3


1)



2) $+45^\circ$ tek5.csv tek5.png



-45° tek6.csv tek6.png



5.3

5.3.1

- 1) The Transmitted Signal is only visible when we around the point of V_{pp} meaning we are at the linear area of the Intensity-Voltage-Ratio which roughly follows a \sin is periodically
- We saved the picture


0.56 kV tek7.png

0.74 kV tek8.png

13.12.18

6.2 Atom

1) After adjusting the set up to be symmetrical in intensity



we start to measure.


The Amplitude Voltage is $14V \pm 0.1V$ $d = 193.5 \text{ cm} \pm 0.2 \text{ cm}$

controls frequency $\rightarrow V_{pp} [\text{V}]$	$f [\text{MHz}]$	$\Delta x_1 [\text{cm}]$	$\Delta x_2 [\text{cm}]$
9	$(85.0 \pm 0.5) \text{ MHz}$	(5.2 ± 0.1)	(10.0 ± 0.1)
8	(91.0 ± 0.5)	(5.4 ± 0.1)	(10.7 ± 0.1)
7	(97.5 ± 0.3)	(5.8 ± 0.1)	(11.5 ± 0.1)
6	(104.1 ± 0.3)	(6.1 ± 0.1)	(14.8 ± 0.1)
5	(109.0 ± 0.4)	(6.5 ± 0.1)	(13.0 ± 0.1)
4	(115.6 ± 0.4)	(6.8 ± 0.1)	— can't see both second maximum
3	(121.0 ± 0.4)	(7.1 ± 0.1)	—
2	(126.5 ± 0.2)	(7.4 ± 0.1)	—
1	(132.3 ± 0.2)	(7.8 ± 0.1)	—

Measure of intensity of 0-1- maximum				Sound wave Amplitude
$V_{pp} [\text{V}]$	$f [\text{MHz}]$	$I_0 [\text{V}]$	$I_1 [\text{V}]$ (left Maxima)	A
9	(85.2 ± 0.2)	(3.35 ± 0.01)	(1.54 ± 0.01)	(17.5)
8	(91.1 ± 0.2)	(3.62 ± 0.01)	(1.08 ± 0.01)	(22.5)
7	(97.6 ± 0.2)	(4.06 ± 0.01)	(0.51 ± 0.01)	(30.0)
6	(104.2 ± 0.2)	(4.50 ± 0.01)	(0.334 ± 0.003)	(27.9)
5	(109.9 ± 0.1)	(4.84 ± 0.02)	(0.150 ± 0.003)	(24.1)
4	(115.4 ± 0.3)	(5.04 ± 0.02)	(0.165 ± 0.003)	(24.1)
3	(121.0 ± 0.2)	(5.10 ± 0.01)	(0.188 ± 0.003)	(32.1)
2	(126.6 ± 0.3)	(5.06 ± 0.01)	(0.240 ± 0.004)	(53.1)
1	(132.3 ± 0.3)	(5.15 ± 0.01)	(0.173 ± 0.006)	(24.2)

6.2.2

1. By combining two AOMs perpendicular to each other we expect to see diffraction in both x- and y-direction. i.e. a lattice. We get a pattern like the following.



2. We first change the voltage of the first AOM (we lower the V_{pp}) the spots get apart from each other in x direction. Then we change the voltage of the second AOM (also lower the V_{pp}), the spots get apart from each other in y direction.

We disconnected the cable from AOM modulator to the power supply and instead connected it from AOM modulator to function generator.


We disconnected the cable from AOM modulator LEVEL TO power supply and instead connected it to a sine function generator.

1. It blinks.

Changing wave form doesn't change the behavior of diffraction pattern.



2. If we change the wave form to $\sin u$, the spots expand and shrink around the central spot.

Rectangle form: the spots spring between two positions (block)

Sinus: 
ramp 

diffraction pattern




~~block~~ 
~~block~~ 

triangle form: similar to sinus

